

Viewing Images in your VistaView

Please follow the steps below to learn how to view the images within your VistaView.


Note: Images will be saved once you accept an image on the imaging screen.

You can view patient images in multiple areas of the application:

- In the exam review page after imaging a patient
- In the patient's record
- While transferring files to a computer or other application.

How to view images on the exam review page

1. After imaging a patient and accepting the images, select **Done** in the lower-right corner of the window.
2. This will bring you to an exam review page. You will see thumbnails of the images you accepted. Tap on them to enlarge the image.
3. This will bring you to a screen where you can further examine the image by zooming in and adding a red-free filter if necessary:






How to View Images in the Patient Record

1. Start at the main menu and select **Manage Patients**
2. Select **Patient List**
3. Select **View** next to the patient you would like to review and the date it was taken
4. This will bring you to the exam review page listed previously where you can tap on an image to enlarge it.

How to View Images in the Patient Record

1. Start at the main menu and select **Manage Patients**
2. Select **Transfer Patient Files**
3. Select **Transfer** next to the patient you would like to review.
4. It will ask you to enter in your password to decrypt the images. Enter your password and select **Confirm**.
5. Tap on the folder on the images you would like to review.
6. Upon selecting an image to review, it will ask you to open it in a different application. We recommend using “Photos” or “Gallery” to open images. Select **Just Once**.



If you need any additional assistance, please contact customer support at support@volk.com, or call us at 440-942-6161.