



LENS SPECIFICATIONS

BIO LENSES				
Lens Name	FOV	Image Mag	Working Distance	Lens Diameter
Macula Plus® 5.5	36° / 43°	5.50x	80 mm	63.2 mm
14D	36° / 47°	4.30x	75 mm	57.4 mm
15D	36° / 47°	4.11x	72 mm	57.4mm
20D	46° / 60°	3.13x	50 mm	55.4 mm
Pan Retinal® 2.2	56° / 73°	2.68x	40 mm	57.4 mm
25D	52° / 68°	2.54x	38 mm	50.1 mm
28D	53° / 69°	2.27x	33 mm	45.9 mm
30D Small	44° / 57°	2.09x	31 mm	34.9 mm
30D	58° / 75°	2.15x	30 mm	48.3 mm
40D	69° / 90°	1.67x	20 mm	45.3 mm
Digital Clear Field	55° / 72°	2.79x	37 mm	51.9 mm
Digital Clear Mag	38° / 49°	3.89x	60 mm	51.9 mm
20D ACS®	46° / 60°	3.13x	50 mm	55.4 mm
28D ACS®	53° / 69°	2.27x	33 mm	45.9 mm