



## LENS SPECIFICATIONS

SINGLE-USE				
Lens Name	Product Code	Mirror Angles	Image Mag	Contact Diameter
Capsulotomy	VCAPSD10	N/A	1.57x	14.2 mm
Iridotomy	VIRIDD10	N/A	1.70x	14.2 mm
SLT	VSLTD10	63°	1.00x	12.0 mm
3-Mirror Gonio	V3MIRD10	60°/ 66°/ 76°	1.00x	8.1 mm
4-Mirror Gonio	V4MIRD10	4x63°	1.0x	8.1 mm
20D	V20LCD10	N/A	3.13x	Non-Contact
28D	V28LCD10	N/A	2.27x	Non-Contact