



# LENS SPECIFICATIONS

GONIO LENSES - GLASS (G-SERIES)				
Lens Name	Mirror Angles	Image Mag	Working Distance	Contact Diameter
G-1 Gonio, With Flange	62°	1.50x	Contact	15.0 mm
G-1 Gonio, No Flange	62°	1.50x	Contact	8.4 mm
G-2 Gonio, With Flange	60° / 64°	1.50x	Contact	15.0 mm
G-2 Gonio, No Flange	60° / 64°	1.50x	Contact	8.4 mm
G-3 Gonio, With Flange	60° / 66° / 76°	1.06x	Contact	15.0 mm
G-3 Gonio, No Flange	60° / 66° / 76°	1.03x	Contact	11.4 mm
G-3 Gonio Mini, No Flange	4x64°	1.00x	Contact	9.6 mm
G-4 Gonio, With Flange	4x64°	1.00x	Contact	15.0 mm
G-4 Gonio, No Flange	4x64°	1.00x	Contact	8.4 mm
G-4 High Mag Gonio, With Flange	4x64°	1.50x	Contact	15.0 mm
G-4 High Mag Gonio, No Flange	4x64°	1.50x	Contact	8.4 mm
G-6 Gonio, No Flange	6x63°	1.00x	Contact	8.4 mm

## GONIO LENSES - ACRYLIC

Lens Name	Mirror Angles	Image Mag	Working Distance	Contact Diameter
3-Mirror, No Flange Coated and Uncoated	60° / 66° / 76°	0.90x	Contact	15.3 mm
3-Mirror, ANF+ Flange Coated and Uncoated	60° / 66° / 76°	0.90x	Contact	18.0 mm
Mini 4-Mirror	4x62°	0.90x	Contact	12.75 mm