



VistaView Technical Specifications

Parameter

- Device Model
- Imaging Modality
- Minimum Pupil Size
- Field of View
- Weight
- Focus Adjustment
- Diopter Correction
- Light Source
- Image Resolution
- Image Format
- Dimensions
- Connector Type
- Device Power Rating
- Battery Capacity
- Software
- Operating System
- Compatible Smartphone
- Switching Power Supply

Description

Volk VistaView

Mydriatic

4 mm

Static: 55° horizontal, 42° vertical

560 grams

Auto-focus, manual focus

-15 to +15 diopter

White light LED

28.4 pixels/degree horizontal,
24.8 pixels/degree vertical

JPEG, DNG, DICOM

88 x 202 x 103 mm

USB-C

5V (DC), 2A

4630 mAh

VistaView Smartphone Application

Android OS

Nokia XR20 (5G Network Compatible)

Input: 100-240V (AC), 50/60 Hz

Output: 5V (DC), 2A