



## LENS SPECIFICATIONS

SLIT LAMP LENSES				
Lens Name	FOV	Image Mag	Working Distance	Lens Diameter
60D	68° / 81°	1.15x	13 mm	34.9 mm
78D	81° / 97°	0.93x	8 mm	34.9 mm
90D	74° / 89°	0.76x	7 mm	25.8 mm
Super 66*	80° / 96°	1.00x	11 mm	34.5 mm
Super Field*	95° / 116°	0.76x	7 mm	30.0 mm
Super VitreoFundus*	103° / 124°	0.57x	4-5 mm	26.7 mm
Super Pupil* XL	103° / 124°	0.45x	4 mm	23.6 mm
Digital Wide Field*	103° / 124°	0.72x	4-5 mm	33.0 mm
Digital High Mag*	57° / 70°	1.30x	13 mm	31.1 mm
Digital 1.0x Imaging	60° / 72°	1.00x	12 mm	34.9 mm