

BRAND	MODEL	MERLIN® CLA MOUNTING HARDWARE	REINVERTER
Alcon	LuxOR	Alcon Mounting Kit (11732)	ROLS∞ auto (11363) ROLS∞ manual (11313) ROLS (VROLS-Z)
Haag-Streit / Möller-Wedel	Allegra 90	DOES NOT MOUNT	ROLS∞ auto (11363) ROLS∞ manual (11313) ROLS (VROLS-Z)
Haag-Streit / Möller-Wedel	Allegra 590	DOES NOT MOUNT	ROLS∞ auto (11363) ROLS∞ manual (11313) ROLS (VROLS-Z)
Haag-Streit / Möller-Wedel	Allegra 900	DOES NOT MOUNT	ROLS∞ auto (11363) ROLS∞ manual (11313) ROLS (VROLS-Z)
Haag-Streit / Möller-Wedel	CooperVision	Moeller Wedel Mounting Kit (11410)	ROLS∞ auto (11363) ROLS∞ manual (11313) ROLS (VROLS-Z)
Haag-Streit / Möller-Wedel	Hi-R 900	Moeller Wedel Hi-R 900 Mounting Kit (12161)	ROLS∞ auto (11363) ROLS∞ manual (11313) ROLS (VROLS-Z)
Haag-Streit / Möller-Wedel	Hi-R 900 SL	DOES NOT MOUNT	ROLS∞ auto (11363) ROLS∞ manual (11313) ROLS (VROLS-Z)
Haag-Streit / Möller-Wedel	Hi-R 900 SLA	DOES NOT MOUNT	ROLS∞ auto (11363) ROLS∞ manual (11313) ROLS (VROLS-Z)
Haag-Streit / Möller-Wedel	EOS 900	DOES NOT MOUNT	ROLS∞ auto (11363) ROLS∞ manual (11313) ROLS (VROLS-Z)
Haag-Streit / Möller-Wedel	Hi-R NEO 900	Moeller Wedel Hi-R 900 Mounting Kit (12161)	ROLS∞ auto (11363) ROLS∞ manual (11313) ROLS (VROLS-Z)
Haag-Streit / Möller-Wedel	Ophthalmic 900	DOES NOT MOUNT	ROLS∞ auto (11363) ROLS∞ manual (11313) ROLS (VROLS-Z)
Haag-Streit / Möller-Wedel	Spectra 900	DOES NOT MOUNT	ROLS∞ auto (11363) ROLS∞ manual (11313) ROLS (VROLS-Z)
Huvitz	HOM-700	Leica Mounting Kit (11409) & Vertical Spacing Kit (11425)	ROLS∞ auto (11364) ROLS∞ manual (11315) ROLS (VROLS-L)
Inami	L-0980A	Inami Mounting Kit (11654)	ROLS∞ auto (11363) ROLS∞ manual (11313) ROLS (VROLS-Z)
Inami	L-0995AZ	DOES NOT MOUNT	ROLS (VROLS-Z)
Leica	M220	DOES NOT MOUNT	ROLS (VROLS-LX)
Leica	M500	Leica Mounting Kit (11409)	ROLS∞ auto (11364) ROLS∞ manual (11315) ROLS (VROLS-L)
Leica	M501	Leica Mounting Kit (11409)	ROLS∞ auto (11364) ROLS∞ manual (11315) ROLS (VROLS-L)
Leica	M520	Leica Mounting Kit (11409)	ROLS∞ auto (11364) ROLS∞ manual (11315) ROLS (VROLS-L)
Leica	M525	Leica Mounting Kit (11409)	ROLS∞ auto (11364) ROLS∞ manual (11315) ROLS (VROLS-L)
Leica	M620	Leica Mounting Kit (11409) & Vertical Spacing Kit (11425)	ROLS∞ auto (11364) ROLS∞ manual (11315) ROLS (VROLS-L)
Leica	M650	Leica Mounting Kit (11409)	ROLS∞ auto (11364) ROLS∞ manual (11315) ROLS (VROLS-L)
Leica	M651	Leica Mounting Kit (11409)	ROLS∞ auto (11364) ROLS∞ manual (11315) ROLS (VROLS-L)
Leica	M655	Leica Mounting Kit (11409)	ROLS∞ auto (11364) ROLS∞ manual (11315) ROLS (VROLS-L)
Leica	M690	Leica M690 Mounting Kit (12159)	ROLS∞ auto (11364) ROLS∞ manual (11315) ROLS (VROLS-L)
Leica	M691	Leica 690 Mounting Kit for CLA (11481)	ROLS∞ auto (11364) ROLS∞ manual (11315) ROLS (VROLS-L)
Leica	M695	Leica 690 Mounting Kit for CLA (11481)	ROLS∞ auto (11364) ROLS∞ manual (11315) ROLS (VROLS-L)
Leica	M820	Leica M820/M840/M841/M844 Mounting Kit (12160)	ROLS∞ auto (11364) ROLS∞ manual (11315) ROLS (VROLS-L)
Leica	M822	Leica M822 Mounting Kit (12171)	ROLS∞ auto (11364) ROLS∞ manual (11315) ROLS (VROLS-L)
Leica	M840	Leica M820/M840/M841/M844 Mounting Kit (12160)	ROLS∞ auto (11364) ROLS∞ manual (11315) ROLS (VROLS-L)
Leica	M841	Leica M820/M840/M841/M844 Mounting Kit (12160)	ROLS∞ auto (11364) ROLS∞ manual (11315) ROLS (VROLS-L)
Leica	M844	Leica M820/M840/M841/M844 Mounting Kit (12160)	ROLS∞ auto (11364) ROLS∞ manual (11315) ROLS (VROLS-L)
Leica	Proveo 8	Leica Mounting Kit (11409) & Vertical Spacing Kit (11425)	ROLS∞ auto (11364) ROLS∞ manual (11315) ROLS (VROLS-L)
Prescott	OMNI-FLEX	Moeller Wedel Mounting Kit (11410)	ROLS∞ auto (11363) ROLS∞ manual (11313) ROLS (VROLS-Z)
Takagi	OM-8	DOES NOT MOUNT	ROLS∞ auto (11363) ROLS∞ manual (11313) ROLS (VROLS-Z)
Takagi	OM-9	Takagi Mounting Kit (11412) & Takagi p/n O06-25	ROLS∞ auto (11363) ROLS∞ manual (11313) ROLS (VROLS-Z)
Takagi	OM-18	Takagi Mounting Kit (11412)	ROLS∞ auto (11363) ROLS∞ manual (11313) ROLS (VROLS-Z)
Topcon	OMS-90	Topcon OMS-90 Mouting Kit (11947)	ROLS∞ auto (11363) ROLS∞ manual (11313) ROLS (VROLS-Z)
Topcon	OMS-600	DOES NOT MOUNT	ROLS∞ auto (11363) ROLS∞ manual (11313) ROLS (VROLS-Z)
Topcon	OMS-610	Topcon Mounting Kit (11413)	ROLS∞ auto (11363) ROLS∞ manual (11313) ROLS (VROLS-Z)
Topcon	OMS-710	Topcon Mounting Kit (11413)	ROLS∞ auto (11363) ROLS∞ manual (11313) ROLS (VROLS-Z)
Topcon	OMS-800 Standard	Topcon Mounting Kit (11413)	ROLS∞ auto (11363) ROLS∞ manual (11313) ROLS (VROLS-Z)
Topcon	OMS-800 Pro	Topcon Mounting Kit (11413)	ROLS∞ auto (11363) ROLS∞ manual (11313) ROLS (VROLS-Z)
Zeiss	OPMI 1	Zeiss Mounting Kit (11408)	ROLS∞ auto (11363) ROLS∞ manual (11313) ROLS (VROLS-Z)
Zeiss	OPMI 6	Zeiss Mounting Kit (11408)	ROLS∞ auto (11363) ROLS∞ manual (11313) ROLS (VROLS-Z)
Zeiss	OPMI 6-CFC	Zeiss Mounting Kit (11408)	ROLS∞ auto (11363) ROLS∞ manual (11313) ROLS (VROLS-Z)
Zeiss	OPMI 6S	Zeiss Mounting Kit (11408)	ROLS∞ auto (11363) ROLS∞ manual (11313) ROLS (VROLS-Z)
Zeiss	OPMI 6-SFC	Zeiss Mounting Kit (11408)	ROLS∞ auto (11363) ROLS∞ manual (11313) ROLS (VROLS-Z)
Zeiss	OPMI 6-SFR	Zeiss Mounting Kit (11408)	ROLS∞ auto (11363) ROLS∞ manual (11313) ROLS (VROLS-Z)
Zeiss	OPMI 6-CH	Zeiss Mounting Kit (11408)	ROLS∞ auto (11363) ROLS∞ manual (11313) ROLS (VROLS-Z)
Zeiss	OPMI CS base w/ Retrolux CS optics	Zeiss Mounting Kit (11408)	ROLS∞ auto (11363) ROLS∞ manual (11313) ROLS (VROLS-Z)
Zeiss	OPMI CS base w/ Retrolux 2 optics	DOES NOT MOUNT	ROLS∞ auto (11363) ROLS∞ manual (11313) ROLS (VROLS-Z)
Zeiss	OPMI CS base w/ Retrolux 3 optics	DOES NOT MOUNT	ROLS∞ auto (11363) ROLS∞ manual (11313) ROLS (VROLS-Z)
Zeiss	OPMI CS base w/ Retroskop optics	DOES NOT MOUNT	ROLS∞ auto (11363) ROLS∞ manual (11313) ROLS (VROLS-Z)
Zeiss	OPMI Lumera	Zeiss Mounting Kit (11408)	ROLS∞ auto (11363) ROLS∞ manual (11313) ROLS (VROLS-Z)
Zeiss	OPMI Lumera i	Zeiss Mounting Kit (11408)	ROLS∞ auto (11363) ROLS∞ manual (11313) ROLS (VROLS-Z)
Zeiss	OPMI Lumera T	Zeiss Mounting Kit (11408)	ROLS∞ auto (11363) ROLS∞ manual (11313) ROLS (VROLS-Z)
Zeiss	OPMI Lumera 300 Basic	DOES NOT MOUNT	ROLS∞ auto (11363) ROLS∞ manual (11313) ROLS (VROLS-Z)
Zeiss	OPMI Lumera 300 Advanced	Zeiss Mounting Kit (11408)	ROLS∞ auto (11363) ROLS∞ manual (11313) ROLS (VROLS-Z)
Zeiss	OPMI Lumera 700	Zeiss Mounting Kit (11408)	ROLS∞ auto (11363) ROLS∞ manual (11313) ROLS (VROLS-Z)
Zeiss	OPMI MD	Zeiss Mounting Kit (11408)	ROLS∞ auto (11363) ROLS∞ manual (11313) ROLS (VROLS-Z)
Zeiss	OPMI MDI	Zeiss Mounting Kit (11408)	ROLS∞ auto (11363) ROLS∞ manual (11313) ROLS (VROLS-Z)
Zeiss	OPMI MDO	ZeissMD Mounting Kit for CLA (11553)	ROLS∞ auto (11363) ROLS∞ manual (11313) ROLS (VROLS-Z)
Zeiss	OPMI MDU	ZeissMD Mounting Kit for CLA (11553)	ROLS∞ auto (11363) ROLS∞ manual (11313) ROLS (VROLS-Z)
Zeiss	OPMI Vario	Zeiss Mounting Kit (11408)	ROLS∞ auto (11363) ROLS∞ manual (11313) ROLS (VROLS-Z)
Zeiss	OPMI VISU 140	DOES NOT MOUNT	ROLS∞ auto (11363) ROLS∞ manual (11313) ROLS (VROLS-Z)
Zeiss	OPMI VISU 150	Zeiss Visu Mounting Kit (12163)	ROLS∞ auto (11363) ROLS∞ manual (11313) ROLS (VROLS-Z)
Zeiss	OPMI VISU 160	Zeiss Visu Mounting Kit (12163)	ROLS∞ auto (11363) ROLS∞ manual (11313) ROLS (VROLS-Z)
Zeiss	OPMI VISU 200	Zeiss Visu Mounting Kit (12163)	ROLS∞ auto (11363) ROLS∞ manual (11313) ROLS (VROLS-Z)
Zeiss	OPMI VISU 210	Zeiss Visu Mounting Kit (12163)	ROLS∞ auto (11363) ROLS∞ manual (11313) ROLS (VROLS-Z)