

# TECHNICAL SPECIFICATIONS

<b>Camera Sensor Resolution</b>	5 Mega Pixel
<b>Optical Resolution</b>	> 60 linear pixels/mm (ISO10940:2009)
<b>Focus Range</b>	Autofocus, -20 to +20D
<b>Minimum Pupil Size</b>	3 mm
<b>Field of View</b>	50 degrees (total FOV 50x40)
<b>Display</b>	4" TFT-LCD, 800x480
<b>Memory</b>	8 GB   4,000+ images
<b>Battery</b>	Rechargeable 3.63V/2600mAh Li-Ion Battery
<b>Image and Video Format</b>	JPEG, DICOM, MPEG-4/MPEG-1
<b>Connectivity</b>	WLAN, USB
<b>Type of Photography</b>	Color, Red-Free, Infra Red, Low-Red
<b>Dimensions</b>	<ul style="list-style-type: none"><li>- Handset - 122(w) x 202(h) x 98 (l) mm, 526g</li><li>- Retina module - 68 (w) x 74 (h) x 160 (l) mm, 327g</li><li>- Anterior module - 71 (w) x 80 (h) x 79 (l) mm, 105g</li></ul>



PICTOR Prestige™

